

LymphaTouch model comparison

Features	LymphaTouch Pro (LT02)	LymphaTouch (LT01)
Maximum pressure	400mmHg	250mmHg*
Treatment programs built in	Lymphedema and Pre/Post-op programs	None *
Materials	Industrial grade plastic handle, steel frame, thicker cable. Durable materials support professional and frequent usage.	Plastic handle and frame for lighter casual use
Handle design	Detachable handle for easier transport and maintenance. Supports optimal grip for conducting manual therapy on customers	Integrated handle. Design supports self-application in home use.
Battery powered	Yes	Yes
Weight	3.1 kg / 6.8 lbs	1.1 kg / 2.4 lbs

*Note: Sport Upgrade and built in programs of the LT01 have been discontinued