

NEW PRODUCT



LYMPHATOUCH®



LYMPHATOUCH PRO

With its combination of negative pressure and mechanical vibration, the LymphaTouch Pro device proves highly effective in softening fibrotic and scar tissue adhesions. This capability enhances tissue elasticity, flexibility, and joint range of motion. As an excellent choice for muscle maintenance and massage treatments, it has the capability to deliver up to 400mHg of negative pressure. It imparts a high-intensity mobilization effect on tissues, thereby enhancing the range and effectiveness of manual therapy.



Key Features:

- Efficient medical device that elevates manual therapy to a new level
- Combines negative pressure and mechanical high-frequency vibration
- Activates the lymphatic circulation and soft tissues in the treated area
- Addresses various acute injuries and musculoskeletal conditions



SKU	Item Name	Specifications
7028016	LymphaTouch Pro Therapy	LT02



Scan to view product online