



## DEQ THUMB ORTHOSIS UNIVERSAL

A versatile solution for thumb-related issues, featuring trimmable straps for customized fit and bendable stays for adaptability. Tailor the fit to each patient's unique needs, ensuring a secure and personalized experience. The open-stay reduces pressure over the first dorsal compartment, alleviating discomfort associated with thumb conditions.

### Key Features:

- Low profile, functional, comfortable
- Open-stay reduces pressure, preventing pain
- Trimmable straps and bendable stay allow custom fit to patient
- Universal size reduces inventory needs



SKU	Item Name	Specifications
7028261	One Size, Univ.	Wrist Circ. 4.5"-9.5"



Scan to view  
product online