




High-Definition Wireless Ultrasound Imaging

You Can Trust for Your Specialty

Meet the world's only 3rd generation wireless ultrasound. Now 30% smaller and lighter. Clarius HD3 delivers best-in-class portable ultrasound for your specialty, with an easy-to-use app powered by artificial intelligence and connected to the cloud.



Improve Diagnostic Confidence and Procedural Safety



Clarius C3 HD3
Curvilinear Scanner

Frequency: 2-6 MHz
Max Depth: 40 cm
Applications: Abdominal, Lung, OB/GYN



Clarius C7 HD3
Microconvex Scanner

Frequency: 3-10 MHz
Max Depth: 18 cm
Applications: Abdominal, Pediatrics, Small Parts, Speech Therapy



Clarius PA HD3
Phased Array Scanner

Frequency: 1-5 MHz
Max Depth: 40 cm
Applications: Cardiac, Lung, Vascular Access



Clarius L20 HD3
Ultra High-Freq. Scanner

Frequency: 8-20 MHz
Max Depth: 4 cm
Applications: Anesthesia, Dermatology, Ophthalmology, Podiatry, Rheumatology, Small Parts



Clarius L15 HD3
High-Freq. Linear Scanner

Frequency: 5-15 MHz
Max Depth: 7 cm
Applications: Breast, MSK, Nerve, Plastic Surgery, Rheumatology, Small Parts, Superficial, Thyroid




Clarius L7 HD3
Linear Scanner

Frequency: 4-13 MHz
Max Depth: 11 cm
Applications: Lung, MSK, Ocular, Vascular




Clarius EC7 HD3
Endocavity Scanner


Frequency: 3-10 MHz
Max Depth: 15 cm
Applications: Early OB, IVF, Pelvic, Urology



Abdomen (C3)




Sag Bladder Peds (C7)



Cardiac (PA)




Early OB (EC7)



Facial Artery and Vein (L20)



Patellar Tendon (L15)



PEC Block (L7)

High-Resolution Ultrasound in a Portable Wireless Scanner



High-Definition

Octal Beamforming technology and 192 piezoelectric elements deliver high performance and resolution.



Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.



Easy To Use

Friendly AI-powered app for your Apple or Android that automatically delivers the best imaging.



Specialized

AI-powered with specialized presets, automated workflows, and advanced specialty packages.



Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!




Clarius Cloud and DICOM

Save images to DICOM or manage your exams anywhere, anytime with unlimited Clarius Cloud storage.




Imaging Modes


Take full advantage of our 192 piezoelectric crystals and octal beamforming technology to power imaging modes never before seen in handheld ultrasound. Get the clinical insights you need without the complexity of a large system.




Pulsed Wave Doppler



Needle Enhance



Elastography



Advanced Software

Refine and personalize your experience with Advanced Software Packages powered by artificial intelligence and purpose-designed for your specialty. From dedicated presets to comprehensive measurements, these packages fine tune imaging and streamline workflows.

 **Aesthetics Package**



 **MSK Package**

 **Primary / Critical Care Package**

 **Breast Package**

 **Obstetrics Package**

 **Vascular Package**



Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than three million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

Clarius Mobile Health

2980 Virtual Way #205,
Vancouver, BC, V5M 4X3
Canada
hello@clarius.com
Phone 1-778-800-9975
www.clarius.com