

# Product Specification Sheet



## Specifications:

Vendor Parts Number: M30MP

### Construction:

- 1:1 weight-to-resistance ratio weights stack.
- 30 lb. weight stack option.
- 1.25 lbs. starting resistance.
- 2" interval rope adjustment along a 67.5" rope adjustment bar.
- TIG welded & powder coated pulley frame.
- Magnetic pin select weight stack.
- Standard accessory: 1x pulley handle.
- Mobility: Tilt & Go feature.
- Assembly: None required.

### Product Dimensions:

- M30MP: 26"(W)x 32"(D)x 68"(T), 75 lbs.

### Shipping Dimensions:

- M30MP: 30"(W)x 36"(D)x 72"(T), 125 lbs. (Class 250)

SHIPS VIA TRUCK 

### Pulley Accessories:

- Straps (Items # 2020, 2020v, 2021, 2022, 2024IM & 2025).
- Pulley Handle (Item # R282H834).
- Bars (Items # 13107 & 13108).
- Support Bar (Item # VH8083).
- Gantry (Item # 15106).
- Eccentric Pulley (Item # SE11579).
- Negative Weight Set (Item # SE9096).
- 3x 0.5 lb. Add-On Weights (Item # 3111).

## Functions:

- Offers light resistance to initiate soft tissue repair and pain reduction during the acute to subacute phase.
- Offers active therapy ranging from active assistive to active resistive in the available range of motion.
- It's size, light weight and "tilt & go" function makes it ideal for active treatment in private rooms and gyms.
- Numerous accessories are available to adjust the exercise dosage to tolerance and individual needs



**Tilt & Go Feature**

RehabPro, Inc.

Phone: 602-317-7207.

E-mail: [sales@rehabpropulleys.com](mailto:sales@rehabpropulleys.com)

Designed by Physical Therapists. Manufactured in the USA.