

SAUNDERS CERVICAL TRACTION

Safe, Effective Home Treatment

SAFE

A proven alternative to surgery for some patients

- Traction forces are directed to the back of the head to prevent compression of the TMJ
- Easy-to-read pressure gauge provides total control of the force at all times
- Blow-off valve limits the amount of force that can be applied

EFFECTIVE

A cost-effective option to continuous clinical care

- Innovative design replicated clinical traction and promotes consistency of treatment
- Positioning alternatives allow varying treatment options for the clinician
- Patented pneumatic system ensures reliable force accuracy and leak-free performance

EASY TO USE

A portable treatment solution that increases patient compliance and satisfaction

- Patented self-adjusting neck wedges ensure a perfect fit
- Patented pump design incorporates a simple “click in place” mechanism
- Quick and simple patient set-up enables home traction without assistance



Performance Health Item #'s:

081206416 SAUNDERS CERVICAL
TRACTION HOMETRAC DELUXE

921145

CONTACT YOUR SALES REPRESENTATIVE

Sales Rep:

Phone:

Email:

PH Preferred

Visit us at
performancehealth.com/saunders-cervical