

BATH & TOILETING





BATHTUB SECURITY POLE AND CURVE GRAB BAR

The Strader Bathtub Security Pole offers the reassurance needed to stand and move confidently, eliminating the fear of falling or losing balance while getting in and out of the bathtub. Its robust construction and ergonomic design make it an essential companion for anyone seeking an added layer of security in the bathroom. With the Curve Grab Bar and user-friendly features, this tool is committed to prevent accidents in one of the most accident-prone areas in the home.

Key Features:

- 3 products-in-1 – grab bar for transferring in and out of the shower, used for stability while showering, and can be rotated to assist in getting on and off the toilet
- Clamps to 3"-7" edge of most bathtubs, and tension mounts to the ceiling without screwing or drilling
- Fits 7' – 9' ceilings
- 360° Rotating grab bar locks into



SKU	Item Name	Specifications
7027485	Bathtub Security Pole and Curve Grab Bar	Fits Bathtub Heights: 13-21 inches



Scan to view product online