

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

METRON2

WOOD TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 550 lbs.
- 30"W x 78"L x 30"H
- H-brace construction
- Maple hardwood legs with hardwood frame
- Heavy duty vinyl with bacterial protection
- 2" High density foam cushion
- Rounded corners
- Edge banding
- Some assembly required

TABLE CONFIGURATIONS:

- FMV-9235**
Flat Cushion
30"W x 78"L x 30"H

Qty _____

- TMV-9235**
Lift Back Cushion
30"W x 78"L x 30"H

Qty _____



FMV-9235



TMV-9235

OPTIONS:

- | |
|---|
| <input type="radio"/> MV-BottomShelf - Bottom Shelf |
| <input type="radio"/> MV-Cabinet - Cabinet |
| <input type="radio"/> MV-NoseHole - Nose Hole |
| <input type="radio"/> MV-PaperCutter - Paper Roll Holder (On Foot End with Paper Cutter on Head End) |

STANDARD UPHOLSTERY CHOICES:

- | | |
|---------------------------------------|---------------------------------------|
| <input type="radio"/> Royal Blue | <input type="radio"/> American Beauty |
| <input type="radio"/> Imperial Blue | <input type="radio"/> Lark |
| <input type="radio"/> Blue Ridge | <input type="radio"/> Burgundy |
| <input type="radio"/> Deep Violet | <input type="radio"/> Dove |
| <input type="radio"/> Forest | <input type="radio"/> Charcoal |
| <input type="radio"/> Mandarin Orange | <input type="radio"/> Mocha |
| <input type="radio"/> Black | <input type="radio"/> Putty |

metron™ **VALUE**

a Performance Health Brand

Toll Free: 800.323.5547 | Fax: 800.547.4333
www.performancehealth.com