

AIDS TO DAILY LIVING




RXCRUSH PILL CRUSHING SYSTEM

RxCrush is designed by a nurse for a nurse. A pill crushing system that delivers pulverized and dissolved medication to a feeding tube. Pills are first crushed and then finely pulverized by the roller within the pill pouch. The resulting powder is ready to be dissolved in liquid.

Key Features:

- Pulverizes and dissolves medication
- Connects to ENFit®
- 2-stage crushing and rolling
- Zip-closed pouches



Scan to view
product online

SKU	Item Name	Specifications
7035453	Pill Crusher II	Pill Crushing System
7035454	Pill Pouch, ENFit®	3"x5", 50ml
7035455	Pill Pouch, ENFit®	2"x3", 30ml
7035456	Sealed Tear-off Pill Pouch	2"x5", 30ml
7035457	Sealed Tear-off Pill Pouch	3"x7", 50ml
7035458	Open Pill Pouch	2"x5"
7035459	Female ENFit® Pill Pouch Caps	ABS Plastic
7035460	Female ENFit® Pill Pouch Caps	Vinyl