



THERMOSKIN ELASTIC ANKLE WRAP

This Elastic Ankle Wrap by Thermoskin is a woven sock wrap support that fits snugly to the ankle. The adjustable hook and loop strap allows additional compression either medially or laterally where needed. The design of this product provides a cool to wear build and allows free ventilation. For sizing, measure around the foot under the arch.

Key Features:

- Used for compression and pain alleviation purposes
- Available in two size options:
- Hand wash in warm water using mild soap and air dry for best results
- Designed for multiple uses



SKU	Item Name	Specifications
081045871	Elastic Ankle Wrap	Regular, 6.75"-9.5"
081045889	Elastic Ankle Wrap	XL, 9.5"-12.25"



Scan to view
product online