

Ladies' Sheer Moderate Support Knee High Style # 1651 (15 – 20 mmHg)

Description: Ladies' sheer moderate support knee high. Graduated compression. Physician recommended. Relieves tired, aching and swollen legs, often symptoms of varicose veins.

Fiber Content: Nylon 80%
Spandex 20%

Compression in mmHg: Ankle 16 mmHg
(Target) Calf 11 mmHg

- Features:**
- Fashionable sheer appearance for work or leisure
 - Two way stretch design provides extra comfort
 - Reciprocated heel and toe
 - Wide comfort top welt
 - Comfortable non-binding top
 - Graduated compression helps shape-up legs
 - Latex free

Sizes:		Small	Medium	Large	X-Large
Colors:	Beige	0815-68-682	0815-68-690	0815-68-708	0815-68-716

Size	Ankle	Calf	Calf Length
SMALL	7" - 8"	10 1/2" - 14 1/2"	TO 16"
MEDIUM	8" - 9 1/2"	11 1/2" - 15 1/2"	TO 17"
LARGE	9 1/2" - 11"	12 1/2" - 17"	TO 18"
X-LARGE	11" - 12 1/2"	13 1/2" - 18"	TO 19"





**PATTERSON
MEDICAL**

www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada call 1-877-566-8611