

Ladies' Sheer Firm Support Knee High Style # 1670 (20 – 30 mmHg)

Description: Ladies' sheer firm support knee high. Graduated compression. Physician recommended. Relieves tired, aching and swollen legs, mild varicosities and edema. Effective as post sclerotherapy treatment.


Fiber Content: Nylon 80%
Spandex 20%

Compression in mmHg: Ankle 23 mmHg
(Target) Calf 16 mmHg

- Features:**
- Fashionable sheer appearance for work or leisure
 - Two way stretch design provides extra comfort
 - Reciprocated heel and toe
 - Comfortable non-binding top
 - Graduated compression helps shape-up legs

Size	Small	Medium	Large	X-Large	2X-Large	3X-Large
Beige	0815-68-724	0815-68-732	0815-68-740	0815-68-757	0815-68-823	0815-68-831

Size	Ankle	Calf	Length
Small	7-8"	10 1/2" – 14 1/2"	To 16"
Medium	8-9 1/2"	11 1/2" – 15 1/2"	To 17"
Large	9 1/2 - 11"	12 1/2" – 17"	To 18"
X-Large	11-12 1/2"	13 1/2" – 18"	To 19"





**PATTERSON
MEDICAL**

www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada call 1-877-566-8611