



SAMMONS PRESTON UNIVERSAL THUMB SPICA SPLINT

The design of the Sammons Preston Universal Thumb Spica Splint stabilizes the thumb and utilizes soft, breathable material to improve patient comfort. An adjustable stay increases support along the thumb but doesn't inhibit movement. The fingers and palm maintain a full range of motion because of the low-profile design. This brace offers a great balance of support and freedom.

Key Features:

- Stabilizes and supports the thumb after injury or strain
- A low-profile design keeps the brace from limiting motion or causing discomfort
- Choose between a short or long version, depending on the type of support required
- The long version is 9". The short version is 7.5"



SKU	Item Name	Specifications
081607423	Short	7.5" Long
081607431	Long	9" Long



Scan to view
product online