

BAILEY

EXTREME

1000 LB Treatment Tables

Model 4510 Bariatric Treatment Table Features:

- 34"W x 28"H x 78"L; optional sizes available.
- 2" thick top, vinyl upholstery in US 387 Camel, standard.
- Choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes.
- Standard H-Brace; *Optional Wide Plain Shelf available as Model 4511.*
- **Six legs and special frame design allows 1000 lb patient capacity!**
- Ships Fully Assembled!
- Ship weight 200 lbs.



Model 4511 Bariatric Treatment Table
[with Wide Plain Shelf]

Model 4515 Adjustable Back Bariatric Treatment Table Features:

- Same features as Plain Table, with added Adjustable Back
- Adjustable Back features positive locking in all positions.
- Standard H-Brace, *Optional Wide Plain Shelf available as Model 4516.*
- **Six legs and special frame design allows 1000 lb patient capacity!**
- Ships Fully Assembled!
- Ship weight 200 lbs.

Model 4515 Adjustable Back Bariatric Treatment Table
[Shown in EX 269 Garnet]



See our full catalog of products at
www.baileymfg.com