



ROLYAN WRIST SUPPORT

The Rolyan Wrist Support is designed with a conforming, creep-resistant elastic that prevents bunching, ensuring a secure and comfortable fit throughout wear. Whether for rehabilitation or daily use, this wrist support combines durability with a tailored fit to help stabilize the wrist without sacrificing comfort or flexibility.

Key Features:

- Conforming, creep-resistant elastic
- Convenient hook & loop closure
- Removable, Adjustable Aluminum Stay
- Available in multiple sizes to ensure a comfortable, customized fit



Scan to view
product online

SKU	Item Name	Specifications
081285345	ROLYAN WRIST SUPPORT RIGHT SM	5.5"-6.5"
081285352	ROLYAN WRIST SUPPORT LEFT SM	5.5"-6.5"
081285360	ROLYAN WRIST SUPPORT RIGHT MED	6.5"-7.25"
081285378	ROLYAN WRIST SUPPORT LEFT MED	6.5"-7.25"
081285386	ROLYAN WRIST SUPPORT RIGHT LG	7.25"-8"
081285394	ROLYAN WRIST SUPPORT LEFT LG	7.25"-8"
081285402	ROLYAN WRIST SUPPORT RIGHT XL	8"-8.75"