

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

ESSENTIAL THERA-P™ BARIATRIC ELECTRIC TREATMENT TABLE

STANDARD FEATURES:


- Load capacity: 650 lbs.
- 34" or 36"W x 77"L x 17-35"H
- Durable and stylish steel base (EPT9020BN)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 34" or 36" Wide cushions in 4 different configuration options
- 2" High density foam cushions with rounded corners
- Individual locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- ADA Compliant

NOTE: International (240V) option available



EPT3012

TABLE CONFIGURATIONS:



EPT3011
1 Section
34" Width

QTY: _____

EPT3611
1 Section
36" Width

QTY: _____

EPT3012
2 Section
34" Width

QTY: _____

EPT3612
2 Section
36" Width

QTY: _____

EPT3013
3 Section
34" Width

QTY: _____

EPT3613
3 Section
36" Width


QTY: _____

EPT3015
5 Section
34" Width

QTY: _____

EPT3615
5 Section
36" Width

QTY: _____




EPT3011 or EPT3611



EPT3012 or EPT3612



EPT3013 or EPT3613



EPT3015 or EPT3615



Toll Free: 800.743.7738 | Fax: 605.882.8398
phsmedicalsolutions.com

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

OPTIONS:

OPTIONS:

- EPT9020BN** - 650 lb Lifting Capacity
- S1077** - Dual Foot Control Pedal
- S1078** - Dual Foot Bar Pedal
- PT10001** - Nose Hole Plug
- S1154** - (1) Battery and (1) Wall Charger
- S1155** - (2) Batteries and (1) Wall Charger



Dual Foot Control Pedal (x2) (S1077)



Dual Foot Bar Control (S1078)

STANDARD FRAME COLORS:

-  White

CUSTOM BASE COLORS:

- | | |
|---|--|
| <input type="radio"/>  Black | <input type="radio"/>  Green |
| <input type="radio"/>  Light Gray | <input type="radio"/>  Forest Green |
| <input type="radio"/>  Charcoal | <input type="radio"/>  Purple |
| <input type="radio"/>  Bright Blue | <input type="radio"/>  Yellow |
| <input type="radio"/>  Royal Blue | <input type="radio"/>  Orange |
| <input type="radio"/>  Navy Blue | <input type="radio"/>  Brown |
| <input type="radio"/>  Red | <input type="radio"/>  Sky Blue |
| <input type="radio"/>  Burgundy | |

STANDARD UPHOLSTERY CHOICES:

- | | |
|---|---|
| <input type="radio"/>  Royal Blue | <input type="radio"/>  American Beauty |
| <input type="radio"/>  Imperial Blue | <input type="radio"/>  Lark |
| <input type="radio"/>  Blue Ridge | <input type="radio"/>  Burgundy |
| <input type="radio"/>  Deep Violet | <input type="radio"/>  Dove |
| <input type="radio"/>  Forest | <input type="radio"/>  Charcoal |
| <input type="radio"/>  Mandarin Orange | <input type="radio"/>  Mocha |
| <input type="radio"/>  Black | <input type="radio"/>  Putty |

CUSTOM BASE OPTIONS:

- S1178B** - Table Base Graphics
- PT9024-CC** - Pick Custom Base Color
- S2000B** - Custom Table Base Powder Coat Color and Table Base Graphics



Toll Free: 800.743.7738 | Fax: 605.882.8398
phsmedicalsolutions.com