

SB3T Roll-InBuddy Tilt

For roll in showers



The tilt is 0-30 degrees,
infinite positioning.



Product features:

- Weight capacity - 160kg (350lb)
- Tilt angle 0 – 30 degrees, infinite positioning
- Rolls over most standard toilets
- Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Safety mechanisms
- Footrests - removable and swing away, height adjustable
- Corrosion resistant Aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap & chest belts plus calf strap
- Caster Wheels - 125mm (5") x 4 locking
- Commode system
- Neck rest, fully adjustable
- Tool-free assembly
- Full details see page 31-36

Accessory items: Refer to the accessory pages or on our website

Product Dimensions:

- Length including footrests 1100mm (43.3")
- Length without Footrests 840mm (33")
- Width – over arms 570mm (22.5")
- Width – Inside Arms 475mm (18.7") - smaller users see accessory pages
- Width - Inside Main legs 400mm (15.7")
- Back Width 460mm (18")
- Back Rest Height from Cushion 440mm (17")
- Toilet Height Clearance (maximum) 525mm (20.7")
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 450 W x 420mm D (17.7" x 16.5")
- Seat Cushion Height from Floor 510mm to 610mm (Adjustable) (20" – 24" Adjustable)
- Seat Cushion Opening 175mm W x 335mm L (6.8" x 13")

* To clear some toilets, the chair may need to be raised from the lowest setting. The struts, which control the tilting function, may restrict the maximum toilet clearance on some toilet types to 465mm (18.3")