

Item # 52018 (pair)

PattStrapTM by *fabrifoam*

An effective device for the immediate comfort and relief for the pain of **Patella Tendonitis, Osgood Schlatter's disease & Iliotibial band Syndrom(ITBS).**

Application Instructions

Patella Tendonitis & Osgood Schlatter's:

Place foam side against skin. Wrap firmly applying pressure directly to the patella tendon (below kneecap) Affix with the attached Velcro ® tab.

For Iliotibial Band Syndrome (ITBS)

Place the strap approximately 2 to 3 inches above the knee. Pull the strap firmly to achieve secure, comfortable compression. Affix the attached Velcro® tab.

Care Instructions:

Rinse product in warm water after each use. Gently squeeze out excess water. Lay flat, foam side up, to dry. The use of detergent is not required for everyday care, when needed hand wash with mild fragrance free detergent. Rinse thoroughly after washing.

**DO NOT MACHINE WASH OR DRY!!
THIS WILL DAMAGE YOUR PRODUCT!**

Manufactured by: Fabrifoam Products
900 Springdale Dr. Exton, PA 19341. www.fabrifoam.com

US Patent 5,036838

Item # 52018 (pair)

PattStrapTM by *fabrifoam*

An effective device for the immediate comfort and relief for the pain of **Patella Tendonitis, Osgood Schlatter's disease & Iliotibial band Syndrom(ITBS).**

Application Instructions

Patella Tendonitis & Osgood Schlatter's:

Place foam side against skin. Wrap firmly applying pressure directly to the patella tendon (below kneecap) Affix with the attached Velcro ® tab.

For Iliotibial Band Syndrome (ITBS)

Place the strap approximately 2 to 3 inches above the knee. Pull the strap firmly to achieve secure, comfortable compression. Affix the attached Velcro® tab.

Care Instructions:

Rinse product in warm water after each use. Gently squeeze out excess water. Lay flat, foam side up, to dry. The use of detergent is not required for everyday care, when needed hand wash with mild fragrance free detergent. Rinse thoroughly after washing.

**DO NOT MACHINE WASH OR DRY!!
THIS WILL DAMAGE YOUR PRODUCT!**

Manufactured by: Fabrifoam Products
900 Springdale Dr. Exton, PA 19341. www.fabrifoam.com

US Patent 5,036838