



ReadyWrap Arm			SIZING (in Centimeters)		
Item Code	Item Desc.	Lngh.	Wrist	Elbow	Axilla
RW-UE-CG-S43	Short Small	43	13.5 - 17.0	19.0 - 26.0	22.0 - 29.0
RW-UE-CG-M43	Short Medium	43	15.5 - 20.0	24.0 - 33.0	26.0 - 36.0
RW-UE-CG-L43	Short Large	43	17.5 - 23.0	29.0 - 40.0	31.0 - 43.0
RW-UE-CG-XL43	Short X-Large	43	17.5 - 23.0	30.0 - 43.0	36.0 - 50.0
RW-UE-CG-S46	Average Small	46	13.5 - 17.0	19.0 - 26.0	22.0 - 29.0
RW-UE-CG-M46	Average Medium	46	15.5 - 20.0	24.0 - 33.0	26.0 - 36.0
RW-UE-CG-L46	Average Large	46	17.5 - 23.0	29.0 - 40.0	31.0 - 43.0
RW-UE-CG-XL46	Average X-Large	46	17.5 - 23.0	30.0 - 43.0	36.0 - 50.0
RW-UE-CG-S49	Tall Small	49	13.5 - 17.0	19.0 - 26.0	22.0 - 29.0
RW-UE-CG-M49	Tall Medium	49	15.5 - 20.0	24.0 - 33.0	26.0 - 36.0
RW-UE-CG-L49	Tall Large	49	17.5 - 23.0	29.0 - 40.0	31.0 - 43.0
RW-UE-CG-XL49	Tall X-Large	49	17.5 - 23.0	30.0 - 43.0	36.0 - 50.0