

SammonsPreston®
a Patterson Brand



Two-Tiered Horizontal Bolt Board

Item Sku: 5329

Item Description: High-quality lacquered wood with felt-covered bottom to prevent surface mar. Bolts progress from 1/4" to 1/2", with assorted nuts, bolts and washers. Bolt threads range from course to fine. Board is 11-1/2"L x 7-3/4"W x 4"H. Latex free.

Indications:

The Two-Tiered Horizontal Bolt Board can be used as an evaluation tool as well as therapeutic exercise to increase your clients' ability to perform day-to-day routines and independent living. Helps clients exercise fine motor skills and gains strength in performing common tasks. Great for work simulation tasks. Easy to use in the clinic or for home visits.

Care and Cleaning:

Two-Tiered Horizontal Bolt Board can be wiped clean with a mild detergent – no bleach, no heavy duty cleaning products - to clean. Do not submerge in water. Dry with soft cloth.

Procedures:

To use:

- Start with smallest bolts and screw on nuts and washers to increase range of motion and restore hand, wrist, and forearm use
- Gradually progress to larger bolts to facilitate adapted task performance



Distributed by:

Patterson Medical
28100 Torch Parkway Ste. 700
Warrenville, IL 60555, USA
Phone: 800-323-5547
Direct Phone: 630-393-6000
Fax: 800-547-4333
Direct Fax: 630-393-7600

Patterson Medical
W68 N158 Evergreen Blvd
Cedarburg, WI 53012, USA
Phone: 800-323-5547
Fax: 800-547-4333

Patterson Medical Canada
6675 Mill Creek Drive, Unit 3
Mississauga, ON L5N 5M4
Phone: 800-665-9200
Direct Phone: 905-858-6000
Direct Fax: 905-858-6001

Authorized European Representative
Patterson Medical, LTD.
Nunn Brook Road
Huthwaite, Sutton-in-Ashfield
Nottinghamshire, NG17 2HU, UK
International
Phone: +44 (0) 1623 448 706
Fax: +44 (0) 1623 448 784