

Visit us online
at Skil-Care.com

Skil-CareTM
CORPORATION

29 Wells Avenue
Yonkers, NY 10701

914-963-2040
800-431-2972

SKIN GUARD LEG PROTECTORS

<u>Reorder #</u>	<u>Unit</u>
703073	Pair

Purpose


These soft, durable neoprene pads are designed to protect legs against skin tears and bruises caused by contact with wheelchair footrest support bars.

Application

Wrap pads around wheelchair footrest support bars. Fabric side faces out. Wrap Velcro[®] tabs around pads and secure as tightly as possible.

Easy Care

The Skin Guard Leg Protectors wipe clean with a damp cloth and mild cleanser. Do not launder this product.



Skil-CareTM
CORPORATION

*Advanced design for nursing and therapy.*TM