



Who could benefit from

The Bedside Valet™

Anybody who requires a little help getting in or out of bed. For the individual who strongly values their **independence** and desires to limit their reliance on others, *The Bedside Valet™* will become a valuable aid in maintaining their lifestyle.

Great for people who suffer from:

- Arthritis
- Multiple Sclerosis
- Chronic Back Pain
- Recovering from surgery
 - Stroke
- Hip pain or replacement
 - Amputees
 - Spinal Cord injury
 - Orthopedic injuries




The Bedside Valet™ is for use on *standard home beds, not medical beds.*

It works by securely clamping onto the angle iron of the bed frame with an 8" bracket and tightened by two four armed handles: **No Tools Required!** Four adjustable feet on the floor provide unsurpassed stability and strength. Carpet cleats are included for carpeted floors.

The Bedside Valet fits the following bed frame sizes.

| Frame Style | Minimum Height | Maximum Height |
|-----------------|----------------|----------------|
| 1 1/2" x 1 1/2" | 6 3/4" | 9" |
| 2 x 1 1/2" | 7 1/4" | 9 1/2" |




The Bedside Valet™

by

AmFab™




**Meeting today's
need for
independent living**



- Promote Self Esteem by enabling independent transfer
- Minimize risk of injury to user or care provider
- Four legged support for maximum stability
- Easy to install: *No Tools Required*
- 8" Mounting Bracket for secure installation
- Four Armed handle for easy tightening
- Solid Steel tubing for lasting durability
- Four easily accessible grasping points for convenient use

- **Independence**
- **Mobility**
- **Injury Prevention**



- Attractive Powder Coat finish
- Allows for normal bedding drape
- Adjustable to fit most beds
- 11" Wide frame for one or two handed use
- No open spaces for potential injury

At some point in life many people find themselves in need of assistance getting in or out of bed. This need may arise from an injury, arthritis, stroke, chronic back pain, or any of a multitude of other ailments. No matter the reason, *The Bedside Valet™* can provide that assistance. With four easily accessible grasping points and a design that assures safety and security, the user can rely on *The Bedside Valet™* for that extra help in performing this critical daily activity.

Independence is critical for the self esteem and vitality of an individual. For the **independent senior** who chooses to stay home or in an Assisted Living environment, *The Bedside Valet™* can improve their quality of life and protect their strongly guarded **independence**.

