


Squeeze with Pressure Dot


Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Over-tighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Remove dot from brace. Wrap brace around wrist and fasten. Do not over-tighten



2. Adjust brace so fastener is not placed on the Ulnar Styloid.



3. With the hand in supination, attach dot directly on top of Ulnar Styloid.



4. Pull straps on dot and fasten.

HELY & WEBER

© 2011 Weber Orthopedic Inc.
DBA Hely & Weber All rights reserved.
1185 East Main St., Santa Paula, California 93060
California: 800-221-5465, National: 800-654-3241
International: 805-525-4244, fax: 805-933-2348
U.S. fax: 800-559-5975, www.hely-weber.com

CE

EC REP

MDSS GmbH, Schiffgraben 41
30175 Hannover, Germany


Rev B 1/11

Squeeze with Pressure Dot


Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Over-tighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Remove dot from brace. Wrap brace around wrist and fasten. Do not over-tighten



2. Adjust brace so fastener is not placed on the Ulnar Styloid.



3. With the hand in supination, attach dot directly on top of Ulnar Styloid.



4. Pull straps on dot and fasten.

HELY & WEBER

© 2011 Weber Orthopedic Inc.
DBA Hely & Weber All rights reserved.
1185 East Main St., Santa Paula, California 93060
California: 800-221-5465, National: 800-654-3241
International: 805-525-4244, fax: 805-933-2348
U.S. fax: 800-559-5975, www.hely-weber.com

CE

EC REP

MDSS GmbH, Schiffgraben 41
30175 Hannover, Germany

Rev B 1/11