


DETECTO.
A Division of Cardinal Scale Mfg. Co.

[Trade Shows](#)
[Print Media](#)
[Credit Application](#)
[Jobs](#)
[Dealers](#)
[Contacts](#)
[About Detecto](#)
[View Cart](#)

Quick Links
[Health Care Products](#)
[Food Service Products](#)

BUILDING A BETTER WEIGH



[Enlarge Image](#)
[Price Quote Request](#)

D1130/D1130K	
Bulletin	

D1130/D1130K Personal Floor Scales

A low profile design, heavy duty construction and a non-skid vinyl mat are just a few features of the Detecto D1130 and D1130K scales - proven performers for budget-minded customers. They feature a white, baked-enamel finish and personal colorful indicators to mark one person's or a whole family's weights. The large diameter dial allows for easy readability even for tall people. Designed for personal use, fitness centers, and health care facilities.

D1130/D1130K Specifications	
Capacities:	D1130 - 300 lb x 1 lb D1130K - 130 kg x 500 g
Overall Dimensions:	12" W x 16.75" D x 2.5" H / 305 mm W x 425 mm D x 64 mm H
Platform Dimensions:	11.25" W x 11.75" D / 286 mm W x 298 mm D
Dial Diameter:	7" / 178 mm
Shipping Weight:	9 lb / 4 kg



HEALTH CARE



FOOD SERVICE

QUALITY WEIGHING PRODUCTS SINCE 1900

(800) 641-2008 detecto@cardet.com 203 E. Daugherty Webb City, MO 64870



PRO-HEALTH SCALES

- Both Pounds/Kilograms and Kilograms-Only Versions are Available.
- Easy-to-Read, Large 7-Inch Dial Registers Weight Quickly.
- Features Adjusting Scroll to Set the Zeropoint and Balance.
- Spacious Platform Accommodates Large Shoe Sizes.
- Low-Profile, Compact Design Requires Little Floor Space.

DETECTO MODEL D350 / D350K PROHEALTH SCALES

The precision movement of this sturdy low-profile scale with circular display guarantees exact results and the design with a lot of white and chrome gives the scale a fresh, timeless look. The powder-coated steel casing is built to last a lifetime. The spacious tread area on the

weighing platform is covered in easy-to-clean, black synthetic material for fast wipe-downs. The modern, elegant scale design will make a positive impression in any room. Durable features include reliable corrosion protection by galvanizing the frame and parts.



model D350

D350 SPECIFICATIONS

CAPACITY	350 lb x 1 lb / 160 kg x 1 kg
DIMENSIONS (OVERALL):	13.25" W x 20.25" D x 4.5" H 336 mm W x 514 mm D x 114 mm H
DIMENSIONS (PLATFORM):	10.6" W x 10.6" D 269 mm W x 269 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	17 lb / 8 kg



model D350K

D350K SPECIFICATIONS

CAPACITY	160 kg x 1 kg
DIMENSIONS (OVERALL):	13.25" W x 20.25" D x 4.5" H 336 mm W x 514 mm D x 114 mm H
DIMENSIONS (PLATFORM):	10.6" W x 10.6" D 269 mm W x 269 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	17 lb / 8 kg

DETECTO MODEL D1130 / D1130K PROHEALTH SCALES

A slim, low-profile design, economical pricing, and non-skid vinyl mat are just a few features of the Detecto D1130 and D1130K scales - proven performers for budget-minded customers. They feature a white, baked-enamel finish and personal, colorful indicators to mark one

person's or a whole family's weights. The large 7-inch-diameter dial allows for easy readability even for tall people. Designed for personal use, fitness centers, and health care facilities. Both pounds/kilograms and kilograms-only versions are available.



model D1130

D1130 SPECIFICATIONS

CAPACITY	300 lb x 1 lb
DIMENSIONS (OVERALL):	12" W x 16.75" D x 2.5" H 305 mm W x 425 mm D x 64 mm H
DIMENSIONS (PLATFORM):	11.25" W x 11.75" D 286 mm W x 298 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	9 lb / 4 kg



model D1130K

D1130K SPECIFICATIONS

CAPACITY	130 kg x 500 g
DIMENSIONS (OVERALL):	12" W x 16.75" D x 2.5" H 305 mm W x 425 mm D x 64 mm H
DIMENSIONS (PLATFORM):	11.25" W x 11.75" D 286 mm W x 298 mm D
DIAL DIAMETER:	7" / 178 mm
SHIPPING WEIGHT:	9 lb / 4 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



A Division of Cardinal Scale Manufacturing Co.

203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By: