

# DEFENDER POST-OP KNEE

## Product Features:

- Simple, easy to use hinge settings.
- Drop lock button- 0, 15, 30, 45 degrees.
- Adjustable struts- 17" - 28" length.
- Quick release buckles for easy on & off.
- Cool style offers more breathability.
- Premium style offers full coverage softgoods providing compression and warmth.

## Indications:

- Post-op/post injury immobilization or ROM control.

Product #	Style	Size	Thigh Circum*
0815-63-253	Cool	Regular	Up to 27"
0815-63-261	Cool	XLarge	27 - 35"
0815-63-279	Premium	Regular	Up to 27"
0815-63-287	Premium	XLarge	27 - 35"

\* measure 6" above mid-patella

  
**PATTERSON**  
MEDICAL

[www.pattersonmedical.com](http://www.pattersonmedical.com)

Call toll free:  
**1-800-323-5547**  
In Canada call 1-800-665-9200

Fax toll free:  
**1-800-547-4333**  
In Canada call 1-877-566-8611

**Patterson Medical**  
1000 Remington Blvd., Suite 210  
Bolingbrook, IL 60440

## EASY TO USE!

### Hinge



Unlock, select flexion/  
extension setting and relock

### Drop Lock



Quick release drop  
lock button

### Struts



Press to adjust for length

### Buckles



Quick release for  
easy on and off