


THE BED STEP

Makes Getting In and Out of Bed Easier



For those individuals the need extra support, an optional side panel is available for The Bed Step

The Bed Step

The Bed Step is design to make getting into and out bed easier and safer. The Bed Step is constructed of durable, splinter free, plywood that is finished with an attractive lacquer coating. The two 4 “ high steps are 24” wide and 12” deep and are covered with foot friendly carpet that provides sure footing. The side panel can be placed on either side on the Bed Step (requires assembly) and provides extra stability. An optional side panel kit is available for individuals that need or want extra support when using the Bed Step. The Bed Step is ideal for individuals that have mobility and balance problems. Weight capacity 400 lbs

- Makes getting in and out of bed easier
- Can be place on either side on the bed
- Steps are covered with foot friendly carpet
- Optional side panel kit is available

Item # 0815-70-142 - The Bed Step
Item # 0815-70-159 - Optional Side Panel Kit

To Order Call: 1.800.323.5547

1000 Remington Blvd., Suite 210, Bolingbrook, IL 60440