



### Advantages

- ✓ **EMR/EHR Connectivity Featuring USB Port**  
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB
- ✓ **Hospital Grade Construction**  
Durable steel structure with non-slip platform enables patients to be weighed safely, accurately and comfortably
- ✓ **“Live” Handrails and 180° Swivel Display**  
Handrails may be used during weighing for patient stability and swivel display allows for patient privacy and ease of use
- ✓ **Value-Added Features**  
Operating functions: BMI, Tare, Zero Out, Hold, Live Weight



[www.pattersonmedical.com](http://www.pattersonmedical.com)

Call toll free:  
**1-800-323-5547**  
In Canada call 1-800-665-9200

Fax toll free:  
**1-800-547-4333**  
In Canada call 1-877-566-8611

1000 Remington Blvd., Suite 210 Bolingbrook, IL 60440

*See reverse side for additional specifications*

180° Swivel Display



"LIVE" Handrails  
for patient stability

Master Pack	Shipping Dims. (L x W x H) / Cube	Shipping Wt. (LB)
1	41.89" x 24.65" x 10.43" / 6.23	65 lb



### Specifications

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB
- 3/4" LCD display, with 180° swivel for ease of viewing
- Functions: LB / KG Conversion, LB / KG lockout, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- Handrails may be used during weighing for patient stability
- Platform Size: 22" (w) x 15-3/4" (d)
- (2) Wheels for ease of movement
- 120V adapter included, 6 AA batteries (not included)
- 2 Year Limited Warranty



Reduced or minimum foam in packaging / recyclable cardboard packaging.