



UlcerCARE features and benefits:

- Zipper and pull in outer stocking assists donning
- Provides nominal compression of 40 mmHg* when the medical stocking with the zipper is worn over the compression liner
- Promotes normal lifestyle - regular clothing, bathing and footwear
- Cost-effective choice, JOBST UlcerCARE minimizes nursing care time per patient, per visit
- Easy and convenient to use, JOBST UlcerCARE is designed for maximum patient compliance
- Ready-to-wear therapeutic gradient compression system
- Four convenient sizes available



	With Zipper				
W/O Zipper	Left	Right	Size	Ankle Circumferenc	Calf Circumference
5630-41-01	0815-68-542	0815-68-476	Small	7" - 8 1/4" (17.8-21cm)	11 1/2" - 13 3/4" (29.2-34.9cm)
5630-41-02	0815-68-559	0815-68-484	Medium	8 1/4" - 9 1/2" (21-24.1cm)	13 1/4" - 16" (33.7-40.6cm)
5630-41-03	0815-68-567	0815-68-492	Large	9 1/2" - 10 3/4" (24.1-27.3cm)	15 3/4" - 18 1/2" (40-47cm)
5630-41-04	0815-68-575	0815-68-500	X-L	10 3/4" - 12" (27.3-30.5cm)	18" - 21 1/4" (45.7-54cm)
	0815-68-591	0815-68-518	XX-L	12" - 13" (30.5-33cm)	19" - 22" (48.3-55.9cm)
	0815-68-609	0815-68-526	XXX-L	13" - 14" (33-35.6cm)	19 3/4" - 23" (50.2-58.4cm)
	0815-68-5617	0815-68-534	XXXX-L	14" - 15 1/4" (35.6-38.7cm)	20 1/2" - 23 3/4" (52.1-60.3cm)



www.pattersonmedical.com

Call toll free:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada call 1-877-566-8611