

3pp™ ThumSling NP™



Adjustable stay included
 Apoyo ajustable incluido
 Support réglables inclus
 Formbare Unterstützung inbegriffen

Soft splint for thumb CMC arthritis

Férula suave para la artritis del CMC del pulgar.

Weiche Schiene für Arthritis des CMC des Daumens.

Attelle doux pour l'arthrite de CMC de pouce.



3pp™ ThumSling™ Long



3pp™ Carpal Lift NP™



3pp™ U wrap NP™



3pp™ WristWrap NP™



3pp™ Elbow Wrap™



3pp™ ThumSpica™ Plus

Made in Mexico/ Hecho en México/ Hergestellt in Mexiko/ Fabriqué au Mexique

3pp™ ThumSling NP™

118 Log Canoe Circle,
 Stevensville, MD 21666 U.S.A.


Emergo Australia
 201 Sussex Street
 Darling Park, Tower II, Level 20
 Sydney, NSW 2000 Australia

Emergo Europe
 Molenstrat 15
 2513 BH The Hague
 The Netherlands

Authorized European Representative



3pp™ ThumSling NP™



Size/Talla/ Größe/ Taille

S/M	6"	-	7 1/4"	(15.2 - 18.4 cm)
L/XL	7 1/4"	-	9"	(18.4 - 22.8 cm)

Right/ Derecha/ Rechts/ Droite

S/M P2010-R23
M/L P2010-R34

Left/ Izquierda/ Links/ Gauche

S/M P2010-L23
M/L P2010-L34

To Order/ a comprar/ zum Kauf/ à l'achat

Order Toll Free 888.378.7763 (U.S only)

Local/Int. 410.604.6393

24 hour Fax 410.604.6398

Mon. - Thurs. 8:30 - 5:30 pm, Fri. 8:30 - 5:00 pm EST



www.3pointproducts.com

Material/ Material/ Material/ Matériau:

Nylon, Lycra®, neoprene. No latex.

Nylon, Lycra®, neopreno. No de látex.

Nylon, Lycra®, néoprène. Latex-frei.

Nylon, Lycra®, Neopren. Pas de latex.



Hand or machine wash. Air dry.

Lavar a mano o a máquina. Secar al aire.

Hand oder maschinenwaschbar. Luft trocknen.


Laver à la main ou à la machine. Sècher a l'air.

©2011, 3-Point Products, Inc.

3pp™ ThumSling NP™



Right/ Derecha/ Rechts/ Droite



1



2



3



4



5



6



7



8

www.3pointproducts.com

3-Point Products, Inc


Stevensville, MD U.S.A

888-378-7763

3pp™ ThumSling NP™



Left/ Izquierda/ Links / Gauche



1



2



3



4



5



6



7



8

www.3pointproducts.com

3-Point Products, Inc

Stevensville, MD U.S.A

888-378-7763



3pp™ Removable Stays

3-Point Products is dedicated to providing constant quality improvement. Many of our splints now include a 2" or 3" long flexible stay for additional support wherever it's needed. 3pp stays feature:

- Cushioned fabric on one side and hook closure on the other attach anywhere on soft 3pp splints.
- Flexible stay can be contoured by hand.

3pp™ Splints

Suggested positioning of stays



ThumSling™



ThumSling™ Long



Elbow Wrap™



Wrist Wrap™

Additional stays available in packs of 5.
Call for order information.

888-378-7763

www.3pointproducts.com

©2011 3-Point Products®, Inc. Stevensville, MD 21666 USA