

# Hinged Lateral "J" Stabilizer with Condylar Pads



- J Buttress, in conjunction with Anterior Flap, produces dynamic patella control
- Biaxial Hinges provide medial and lateral stability
- Kuhl™ perforated neoprene allows moisture to escape and the skin to remain dry

J Hinged Brace Sizing:  
Measure at mid-patella with leg fully extended.

XS:	12"-13"	(30.5 cm - 33.0 cm)
Sm:	13"-14"	(33.0 cm - 35.6 cm)
Med:	14"-15"	(35.6 cm - 38.1 cm)
Lg:	15"-17"	(38.1 cm - 43.2 cm)
XL:	16"-18"	(40.6 cm - 45.7 cm)
XXL:	18"-20"	(45.7 cm - 50.8 cm)

Order Number	Description
081578053	KUHL HINGED LATERAL J LEFT XS
081578087	KUHL HINGED LATERAL J LEFT LARGE
081578095	KUHL HINGED LATERAL J LEFT XL
081578103	KUHL HINGED LATERAL J LEFT XXL
081578129	KUHL HINGED LATERAL J RIGHT XS
081578137	KUHL HINGED LATERAL J RIGHT SM
081578145	KUHL HINGED LATERAL J RIGHT MED
081578152	KUHL HINGED LATERAL J RIGHT LARGE
081578160	KUHL HINGED LATERAL J RIGHT XL
081578178	KUHL HINGED LATERAL J RIGHT XXL

Condyle pads provide added comfort and protection



[www.pattersonmedical.com](http://www.pattersonmedical.com)

Call toll free:  
**1-800-323-5547**  
In Canada call 1-800-665-9200

Fax toll free:  
**1-800-547-4333**  
In Canada call 1-877-566-8611

1000 Remington Blvd., Suite 210

Bolingbrook, IL 60440