

Moves on the Sheet...

Not on the Feet!






HEELIFT[®]
Glide
SUSPENSION BOOT



- Suspends the heel in space for a pressure-free environment.
- Soft inner material cradles the calf, redistributing pressure over a wider area and providing interior grip.
- 40% lighter, cooler and less bulky than a pillow-style boot.
- Fabric covering, tricot backing and innovative stitching ensures there are no rough, hard or sharp edges that can damage sensitive skin.



Because proper fit is crucial to the success of your prevention and treatment program, Heelift[®] Glide is available in three sizes.

					
				Smooth	Convuluted
Petite	70 - 130 lbs	6 - 10 in	46 - 64 in	0816-13-900	0816-13-918
Standard	120 - 250 lbs	8 - 14 in	60 - 77 in	0816-13-819	0816-13-801
Bariatric	220 - 600 lbs	12 - 23 in	65 - 80 in	0816-13-892	

What makes the Heelift® Glide the right boot for your patient?



SIMPLE. PROVEN. SOLUTIONS.

Forefoot Strap:
Keeps the boot in place and helps support the forefoot.

Durable Exterior:
Allows the patient's leg to glide easily across bed sheets.

Soft "easy-on, easy-off" Closure Straps:
A simple solution for a custom fit.

Ventilation Holes:
Provide additional air circulation to promote increased patient comfort.

Choice of Two Interiors:
Available with a smooth interior for patients with sensitive skin or edema, or convoluted for better air circulation. (Bariatric available in smooth only)

Heel Suspension Opening:
Allows for visual inspection of the heel.

Friction-free Tricot Backing:
Aids in patient mobility and helps ensure patient safety.

Affixed Elevation Pad:
Lifts the leg to suspend the heel, promoting the ideal environment for an ulcer to heal.

HEELIFT® SAVINGS CALCULATOR

When the dollars make sense...

The cost to treat heel pressure ulcers to discharge can range from \$2,000 (for Stage I) to more than \$70,000 for Stage IV.¹ These costs increase when calculating the cost of treatment to healing.

When used as part of a prevention protocol, Heelift® brand boots can save facilities as much as \$900,000² annually in treatment costs.

Learn how much money your facility could save using Heelift® Glide or Heelift® Suspension Boot at <http://heelift.com/heeliftcalculator.html>

¹ Young ZF, Evans A, Davis J. J Nurs Admin (Jona). Jul/Aug 2003; 33 (7/8): 380-3

² Drennan, Dr. Denis B., "Two Birds With One Stone: How a Professional Offloading Device Can Lower Both Heel Pressure Ulcer Prevalence and Facility Costs," available at http://heelift.com/white_papers/index.html



Reduce costs. Save time. Prevent ulcers.



www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada Fax 1-877-566-8611