



Rolyan® Comfort Collar

Item #081506864 – Small

Item #081506872 - Medium

Item #081506880 - Large

Instructions

1. Wrap collar around neck and fasten securely using hook tab on ends. Once fastened, turn the collar so the closure is at the back.
2. Collar should rest on the clavicle while supporting the neck.
3. To wash, close hook/loop, hand wash or machine wash delicate cycle (in mesh laundry bag preferred) and stretch while wet and let air dry. Do not tumble dry.



WARRANTY

This item is guaranteed to be free from defects in material and/or workmanship under normal use and service for a period of one year from date of purchase. The above guarantee is null and void if unauthorized repair and/or modification is performed on the unit.



Patterson Medical
1000 Remington Blvd.
Suite 210
Bolingbrook, IL 60440
Phone: 800-323-5547
Direct Phone: 630-378-6000
Fax: 800-547-4333
Direct Fax: 630-378-6310



Authorized European Representative
Homecraft Rolyan
A Patterson Medical Company
Nunn Brook Road
Huthwaite
Sutton in Ashfield
Nottinghamshire
NG17 2HU, UK

UK Customer Services
Tel: 08444 124 330
Fax: 08448 730 100

International Customer Services
Tel: +44 1623 448 706
Fax: +44 1623 448 784