


SILIPOS®

Half Gel Coated Body Sleevings



- 10 Foot Rolls – 3”, 4”, 5” Flat Widths
- Cut to Length – Will not Fray
- Coated ½ side with Gel – breathable non site side
- New Beige Color – Lycra/Nylon – will not stretch out
- Good for Hard To reach Sites
- Skin Grafts, Burn, Scar Management, Lymphedema Skin Care

HALF GEL COATED BODY SLEEVING			Circumferencial Measurements	
Product No.	Roll Size	Flat Width	Relaxed	Fully Stretched
08-1451-038	3” x 10’	3”	5.5”	16”
08-1451-046	4” x 10’	4”	7.5”	18”
08-1451-053	5” x 10’	5”	9.5”	20”



Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada fax 1-905-566-9207

www.sammonspreston.com

©2008 Patterson Medical Products, Inc. All rights reserved. Printed in USA.