

TREATMENT/MAT TABLE, BARIATRIC



CLASSIFIED
c UL US IEC60601-1(ed.3)
RM 14971

 750 lb capacity

Motorized Hi-Lo, Scissor

<u>Item No.</u>	<u>Width</u>	<u>Length</u>
TWG8132.BTM	32"	78"
TWG8136.BTM	36"	78"
TWG8142.BTM	42"	78"
TWG8148.BTM	48"	78"

FEATURES:

- 750 lb Weight Capacity, *lifting*
- Motorized Height Adjustable
- 4 Total Lock Casters
- 0 Horizontal Travel Space
- Hand Control
- Steel Top with 2" Removable Mat
- Steel Frame

Height Range: 17 1/2" - 36 3/4"

Finish: Pearl Grey Standard; Porcelain Optional.

Upholstery: Imperial Blue, Regimental Blue, Wedgewood, Blue Ridge, Dove, Graphite, Black, Burgundy, Mauve, Forest, Grotto, Sand, Camel.

Electrical: 120V, 6.5A, 780W, hospital grade plug

Warranty: Ten (10) Year Limited Warranty.

*Nominal Dimensions Stated