



MAT TABLE, BARIATRIC



Motorized Hi-Lo

<u>Item No.</u>	<u>Width</u>	<u>Length</u>
TWG9260.BM	5'	7'
TWG9272.BM	6'	7'

FEATURES:

- 1000 lb Weight Capacity
- 4.0 Safety Factor
- Motorized Height Adjustable
- Motorized Raised Back
- Total Mobility - Roll & Set Leveling Casters; Total Locking
- Handheld Switch
- 2" Upholstered Top - Removable
- Steel Frame with Round Corners
- Bumper Guard Incorporated in Steel Frame
- UL Classified
- Accessories Available: Sit-to-Stand, Bariatric Knee & Chester Stabilizers, Guard Rails, Proprioceptive Activity Components, Mini Cable & Pully Systems, and more.

Height Range:	18" - 38" (<i>nominal</i>)
Raised Back:	Adjustable 0 to 85 degrees
Base:	Heavy Gauge Steel
Finish:	Pearl Grey Standard; Porcelain Optional.
Upholstery:	Imperial Blue, Regimental Blue, Wedgewood, Blue Ridge, Dove, Graphite, Black, Burgundy, Mauve, Forest, Grotto, Sand, Camel.
Electrical:	110V, 13A, 1430W, hospital grade plug.
Warranty:	Ten (10) Year Limited Warranty.

*Nominal Dimensions Stated