



SWEDE-O THERMAL CARPAL TUNNEL GLOVE



Swede-O® Thermal Carpal Tunnel Glove provides comfortable, firm support for weak, injured or aching wrists. Rigid splint helps to prevent wrist movements that can cause Carpal Tunnel Syndrome. Provides heat therapy and even compression to help relieve pain and inflammation.

- Helps prevent wrist movements that can cause Carpal Tunnel Syndrome
- Comfortable, firm support for weak, injured or aching wrists
- Maintains the wrist in a neutral cock-up position
- Anatomically conforming fit and allows for full use of fingers
- Lining captures and retains natural body heat to increase blood flow and aid in the healing process

SIZING GUIDE AND MODEL NUMBERS

Measure the circumference of the hand around the knuckles



LEFT	RIGHT	SIZE	MEASUREMENT	
081694082	081694140	XS	6-6 ³ / ₄ "	15-17 cm
081694090	081694157	S	7-7 ³ / ₄ "	18-20 cm
081694108	081694165	M	8-9"	21-23 cm
081694116	081694173	L	9 ¹ / ₄ -10 ¹ / ₂ "	24-26 cm
081694124	081694181	XL	10 ³ / ₄ -11 ¹ / ₂ "	27-29 cm
081694132	081694199	2XL	11 ³ / ₄ "+	30 cm+

Your Swede-O Thermal Support should be a snug and comfortable fit. Not too tight, not too loose. Use of the correct size is essential for effective results.



Thermal Technology

Swede-O Thermal Carpal Tunnel Glove is designed to provide firm support for weak, injured or aching wrists while reducing the effects of Arthritis and other wrist and hand conditions. The thermal lining captures and retains your natural body heat for effective pain relief and therapy. With unsurpassed quality this thermal support provides comfortable compression and warmth to help ease the pain of sore joints, reduce swelling and associated pain.



www.pattersonmedical.com

Call toll free:
1-800-323-5547

Fax toll free:
1-800-547-4333

Swede-O Thermal Supports are:

