

## **Triaid Cycles for Life – Sizing Guide**

**Terrier** is a rugged tricycle for children from 4 ½ years, this model has 16" wheels and a caliper brake with a parking brake mechanism.

Overall Length 45" (114cm)

Overall Width 29" (74cm)

To suit Inside Leg 16 1/2" (42cm) to 21 1/2" (55cm)

Seat to handlebars 14" (35cm) to 17" (43cm)

Max User Weight 100lbs (45kg)

Wheel Size (Front) 16" x 2.125

Wheel Size (Rear) 16" x 2.125

Seat to ground 19 1/2" (50cm)

**TMX** is a robust tricycle for older children (8 years old), it is fitted with 20" wheels and has a caliper brake with parking brake feature.

Overall Length 59" (150cm)

Overall Width 29" (74cm)

To suit Inside Leg 20" (51cm) to 25" (64cm)

Seat to handlebars 19" (48cm) to 22" (56cm)

Max User Weight 150lbs (68kg)

Wheel Size (Front) 20"

Wheel Size (Rear) 20"

Seat to ground 23 1/2" (60cm)

**Tracker 24"** model is a popular choice for teenagers and adults. It has a robust frame with a low step-through for easy access.

Overall Length 69" (175cm)

Overall Width 32" (81cm)

To suit Inside Leg 25" (64cm) to 30" (76cm)

Seat to handlebars 25" (63cm) to 22" (56cm)

Max User Weight 170lbs (77kg)

Wheel Size (Front) 24" x 2.125

Wheel Size (Rear) 24" x 2.125