

Invacare®  
Rehab  
Shower  
Commode  
Chairs

# Invacare® Rehab Shower Commode Chairs



Yes, you can.™



Yes, you can.™

Invacare's new Rehab Shower Commode Chairs offer aluminum frames and stainless steel hardware that are rust-resistant, making them ideal for use in the shower. These new chairs fit over most standard and elongated bowls for over the toilet use as well. They feature flip-back padded arms that lock down for safety; swing-away front riggings with tool-free adjustable height footrests and a four-position padded seat. The new Invacare® Rehab Shower Commode Chairs offer various seat widths and wheel options for a custom fit along with additional features and benefits.



Model nos. 6895 - 23" treaded urethane tires with 18" seat (shown).  
6795 - 23" treaded urethane tires with 16" seat  
6891 - Four 5" locking casters with 18" seat

**FEATURES**

- Aluminum frame and stainless steel hardware
- Fits over most standard and elongated toilet bowls
- Padded four-position seat
- Flip-back padded arms that lock down for safety
- Swing-away front riggings
- Tool-less adjustable height footrests
- Folds into a compact unit for storage and transport
- Removable pan and pan holder
- Safety strap included
- 23" treaded urethane tire and 5" caster models



Flip-back padded arms that lock down for safety.



Padded four-position seat.



Swing away front riggings with tool-less adjustable height footrests.

**Specifications**

	6895 - 18" seat	6795 - 16" seat	6891 - 18" seat/5" casters
Overall chair width (handrim to handrim)	26 1/2"	24 1/2"	22 1/2" no handrim
Overall back width	18"	16"	18"
Overall chair height (floor to top of push handle)	39"	39"	39"
Overall depth w/o front riggings	32"	32"	29"
Seat width	18 1/4"	16"	18 1/4"
Seat depth	18 1/4"	16"	18 1/4"
Width between arms	18 1/2"	16 1/2"	18 1/2"
Folded width	12"	12"	8"
Seat-to-floor height	22 3/4"	22 1/2"	22 3/4"
Seat-to-armrest height	8"	8"	8"
Back height (top of seat to top of upholstery)	15 1/2"	15"	15 1/2"
Footrest extension min.-max. (from top of seat)	13"-18"	13"-18"	13"-18"
Product weight	41 lb.	39 lb.	32 lb.
Weight capacity	300 lb.	300 lb.	300 lb.
Quantity	1 per carton	1 per carton	1 per carton
Limited warranty	3 years	3 years	3 years

**Clearance for over the toilet use**

Depth (front cross brace to edge of rear tire)	26 1/2"	26 1/2"	23 1/2" no tire
Depth (rear wheel axle bolt to edge of tire)	11 1/2"	11 1/2"	n/a
Height (floor to seat base)	20 1/4"	20 1/2"	20 3/4"
Inside width (between side frames)	16"	14"	16"

**Invacare Corporation**  
www.invacare.com

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare is a registered trademark and 'Yes, you can.', and the blue medallion design are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2004

Invacare Corporation

Form No. 03-101



"Invacare says 'Yes, you can', and I couldn't have said it better myself."