

VELOCITY WITH PATELLA STABILIZATION



US Patent Number 5,759,167

Features and Benefits



- Superior knee brace that provides the control and support your patients need while maintaining the flexibility and the restriction of movement that your rehabilitation protocol requires.
- Constructed in Kuhl™ perforated neoprene which breathes well and offers warmth and maximum compression.
- Range of Motion Hinges feature a polycentric adjustable design to control extension and flexion with comfortable shock absorbing condylar pads and shells which ensure patient comfort and compliance.
- Circumferential straps have elastic hamstring tension relief allowing no interruption in the firing mechanism of the quadriceps muscles while still providing low profile support.
- The brace transforms from left to right by simply moving the patella force strap from the left to right condylar location.
- The brace offers effective support for ligament / meniscus tears while having the patented “u-shaped” patella locator buttress and adjustable patella control of the Shields™ design.

Range-of-Motion Options:

Extension Stops allow 10°, 20°, 30°, 40° and 50° settings, (preset @ 10° upon delivery)

Flexion Stops at 10°, 50°, 60°, 70°, 80°, 90° and 100° increments



Order Number	Description
081578194	VELOCITY PS HINGED KNEE BRACE XS
081578202	VELOCITY PS HINGED KNEE BRACE SM
081578210	VELOCITY PS HINGED KNEE BRACE MED
081578228	VELOCITY PS HINGED KNEE BRACE LG
081578236	VELOCITY PS HINGED KNEE BRACE XL
081578244	VELOCITY PS HINGED KNEE BRACE XXL
081578251	VELOCITY PS HINGED KNEE BRACE 3X

VELOCITY PS SIZING:

Measure 6 inches (15.2cm) above patella.

XS: 5645PS-XS 13"-15" (33.0cm - 38.1cm)
Sm: 5645PS-S 14"-17" (35.5cm - 43.2cm)
Med: 5645PS-M 15"-18" (38.1cm - 45.7cm)
Lg: 5645PS-L 17"-20" (43.2cm - 50.8cm)
XL: 5645PS-XL 19"-22" (48.3cm - 50.8cm)
XXL: 5645PS-XXL 21"-23" (53.3cm - 58.4cm)
XXXL: 5645PS-XXXL 22"-24" (55.8cm - 61.0cm)




PATTERSON
MEDICAL

www.pattersonmedical.com

Call toll free:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada call 1-877-566-8611

1000 Remington Blvd., Suite 210

Bolingbrook, IL 60440