

# MalleoLoc®

Orthosis for stabilization of the ankle.



ProprioPoint  
plantar guide



## Indications


- Early functional treatment for injuries of the lateral malleolar ligaments (and the bifurcate ligaments)
- Capsular ligament strains
- Postoperative rehabilitation
- Chronic ligament insufficiency
- Prophylactic ankle protection

## Features

- Conforms to the outside edge of the foot to counteract supination, even without a shoe.
- Low-profile design, it can easily be worn inside a shoe.
- A special low-fatigue pad makes it extremely comfortable to wear.
- Can be individually shaped.

## Benefits

- The orthosis stabilizes the foot by preventing lateral twisting movements of the ankle without appreciably hindering the normal heel-toe stride.
- MalleoLoc® lies against the antero-lateral aspect of the joint to counteract forward movement of the talus.
- The proprioceptive effects promote active stabilization of the ankle muscles.



Size	1	2	Order number
Heel width in inches	less than 2 1/4	greater than 2 1/4	12013012080+0°
+6 = right, +7 = left, ° = size			