

# ROLYANFIT WRIST BRACE

## Product Features:

- Durable lightweight fabric is lycra lined for breathability and patient comfort.
- Removable aluminum palmar stay provides correct fit and proper support.
- Easy to use D-ring closure.
- 6" and 8" lengths

## Indications:

- Tendonitis, post-cast healing, soft tissue injuries and carpal tunnel syndrome.

6 inch



8 inch



Left	Right	Size	Wrist Circ.
9299-51	9299-47	Small	6¼"-7¼"
9299-52	9299-48	Medium	7¼"-8"
9299-53	9299-49	Large	8"-9"
9299-54	9299-50	Xlarge	9"-10"

Left	Right	Size	Wrist Circ.
9298-87	9298-88	XSmall	5¼"-6¼"
9272-20-02	9272-20-01	Small	6¼"-7¼"
9272-21-02	9272-21-01	Medium	7¼"-8"
9272-22-02	9272-22-01	Large	8"-9"
9272-23-02	9272-23-01	XLarge	9"-10"