





POSITIONING WEDGE

Reorder #	Size	Cover Fabric	Unit
□554010	8"x8"x17"	Nylon With Anti-Slip Mat	Pair
□554020	8"x8"x17"	No Cover	Pair

Purpose

Skil-Care's Positioning Wedge is designed to promote proper hip-spine alignment, reduce sacral pressure to prevent the formation of pressure sores, and aid in the healing of existing sores. Soft and secure, the Position Wedge is more secure than a pillow.



Added Feature The cloth-covered wedges now have anti-slip mats sewn to the bottom. (See illustration) Make certain that the anti-slip mat is positioned so that it is in contact with the bed sheet.

New anti-slip material keeps wedge securely in place.

Easy Care

The nylon cover and anti-slip mat may be cleaned with a damp cloth and mild cleanser. Do not launder cover or foam.