

ROM WALKER

Product Features:

- Provides protected range-of-motion following trauma or post-operative procedures.
- ROM settings between 45° of plantar flexion and 30° dorsi-flexion in 7.5° increments.
- Low profile design provides comfort and stability.
- Shock absorbing insole reduces impact of heel strike.
- Non-skid rocker sole is designed to promote a natural gait.
- Latex free.

Indications:

- Acute ankle sprains, soft tissue injuries, stable fractures, post-operative use following Achilles tendon repair.



Product #	Size	Men's Shoe Size	Women's Shoe Size
5629-60-01	Small	4.5 - 6	5.5 - 7
5629-60-02	Medium	6.5 - 10	7.5 - 11
5626-60-03	Large	10.5 - 12	11.5 - 13