

Comfy™ Grip Hand

Gradual finger extension, contracture management, resting comfort



SKIN SAFE

Anti-microbial material for skin breakdown prevention

COMFORT

Super Soft and removable
Headliner cover is machine washable

ADJUSTABLE

Patented reinforced spine bends to wrist and forearm needs

EXTENSION

Two graduated rolls gradually increase extension of the MCP and PIP joints



PRACTICAL

Holds hand in a functional C position and helps control wrist drop

CONTRACTURE MANAGEMENT

Frame designed to slide between palm and finger of even the most contracted hands

SIDE WINGS

Adjustable wings prevent ulnar and radial deviation

		Wrist to Palmer Creases	Wrist Down Forearm	Splint Length	Patterson #
Left	Large	3 ¼"	7"	10 ¼"	5648-27-03
Right					5648-27-06
Left	Medium	2 ¾"	6 ¼"	9"	5648-27-02
Right					5648-27-05
Left	Small	2"	5"	7"	5648-27-01
Right					5648-27-04