



Deluxe Wrist Drop Orthosis. Lightweight, yet durable Kydex™ plastic can be heat molded for customization and positioning. Holds the hand in extension to support a weak or flaccid wrist. The Universal Cuff, sewn into the palmar strap of the brace, accommodates eating utensils and writing instruments, assisting with daily functions. The removable, washable, foam liner wicks away perspiration and moisture from the skin. Comfortable enough to be worn for extended periods of time. To size, measure width of palm.

Item No.	Side	Size	Price
565879	Right	Large 3½" – 4"	\$68.98
565880	Left	Large 3½" – 4"	\$68.98
565881	Right	X-Large 4" – 4½"	\$68.98
565882	Left	X-Large 4" – 4½"	\$68.98
Item No.	Description		Price
565883	Universal Cuff 90° Angle		\$17.33
565884	Universal Cuff Straight		\$18.81