

Introducing the New ...

# Low PROFILE SI-BELT

**with the Patented Mechanical  
Advantage Pulley System**

The new low profile SI Belt provides stabilization of the sacroiliac (SI) joints through direct circumferential compression utilizing the 3 to 1 ratio mechanical advantage pulley system. The SI belt is prescribed for conservative treatment of sacroiliac strain or arthritis.

**Code L0621 Suggested**



**Cybertech**<sup>™</sup>  
MEDICAL



#### Indications:

- SI joint pain
- Arthritis
- Symphysis pubis dysfunction

#### Features:

- 3:1 mechanical advantage compression power provides unmatched pelvic-sacral support.
- Adjustable sacroiliac orthosis with breathable mesh, comfortable to wear.
- Significantly reduces motion about the sacroiliac joint. Instant pain relief in most cases.
- Customized stabilization pressure specific to the patient.
- Pads are adjustable. Remove both pads for the sleekest fit.

#### Item #:

**0815-56-638** Small  
Fits Hip Circumference Up to 35"

**0815-56-646** Standard  
Fits Hip Circumference 35" – 50"

Excellent for post partum symphysis pubis dysfunction with pads removed



  
**PATTERSON**  
MEDICAL

[www.pattersonmedical.com](http://www.pattersonmedical.com)

Call toll free:  
**1-800-323-5547**  
In Canada call 1-800-665-9200

Fax toll free:  
**1-800-547-4333**  
In Canada call 1-877-566-8611